


# T.O.H.R


Temperature, Oxygen and Heart Rate measuring  
device - Team 2





# INTRODUCTION

T.O.H.R is an advanced health monitoring system designed to measure and track vital physiological parameters in the human body, including heart rate, body temperature, and blood oxygen levels. Inspired by the continuous rhythm of the human heartbeat, this innovative system provides real-time insights into the body's vital signs.



# TABLE OF CONTENTS

## 01

### Objectives

- Measure vital parameters accurately and provide early diagnosis of abnormal conditions.
- Enhance health monitoring, enable timely intervention, and support personalized care for proactive well-being

## 02

### Methodology

- Employing continuous sensing to gather real-time data on vital parameters.
- Utilizing targeted analysis techniques for in-depth examination and insights.

## 03

### Results Analysis

- Visual presentation of analysis results for effective communication with the patient.

## 04


### Conclusions

- Implementation of a soothing routine to promote patient relaxation.
- Continuing the routine until the patient's vital parameters reach a stable state.




# Hardware


»» KY-028



»» DS1307




»» SSD1306




»» MAX32664



»» Buzzer



»» Passive electronics







# Data analysis



# Parameters measured






# Heart rate analysis

Threshold low: 60 bpm  
Threshold high: 100bpm

Temperature measurement is based on the average calculated in 15 seconds of measurements. Correlated by relative accuracy.





# Oxygen analysis

Threshold: 95%



Temperature measurement is based on the average calculated in 15 seconds of measurements.








# Temperature analysis

Threshold: 26°C

Temperature measurement is based on the average calculated in 15 seconds of measurements.



# ANALYSIS RESULTS



## Failure Conditions

- Null value
  - Outliers
  - Confidence
- Null value
  - Outliers
- Environmental temperature

## Elaboration



Each data is associated with an uncertainty value computed on the mean value.  
The relative uncertainty expresses the accuracy of the measure.






# Software protocol





# Finite state automaton





# Finite state automaton - states



# Power consumption analysis

The power consumption of the sensors is based on datasheet informations and it is done in the hypothesis that the sensors are always up and running. This way we can estimate a battery life duration of 13 Hrs 32 Min\*.



\*external battery supplying only sensors with 850mAh capacity at 3.3v



# CONCLUSIONS

Our device is a valid medical aid with good future perspectives as a wearable product that performs continuous monitoring and in depth on demand analysis

